

Micro-spatules



longueur	lame	Ø fil	référence	Prix HT
acier inox AISI 304 18/10				
100 mm	30 x 2 mm	1,0 mm	BC3010	
130 mm	40 x 2 mm	1,0 mm	BC3011	
150 mm	40 x 2 mm	1,0 mm	BC3012	
100 mm	35 x 3 mm	1,5 mm	BC3013	
130 mm	40 x 3 mm	1,5 mm	BC3014	
150 mm	40 x 3 mm	1,5 mm	BC3015	
100 mm	35 x 4 mm	2,0 mm	BC3016	
130 mm	40 x 4 mm	2,0 mm	BC3017	
150 mm	40 x 4 mm	2,0 mm	BC3018	
185 mm	50 x 4 mm	2,0 mm	BC3019	
100 mm	35 x 5 mm	2,0 mm	BC3020	
130 mm	40 x 5 mm	2,0 mm	BC3021	
150 mm	40 x 5 mm	2,0 mm	BC3022	
185 mm	50 x 5 mm	2,0 mm	BC3023	
210 mm	55 x 5 mm	2,0 mm	BC3024	
150 mm	40 x 6 mm	2,0 mm	BC3025	
185 mm	50 x 6 mm	2,0 mm	BC3026	
210 mm	55 x 6 mm	2,0 mm	BC3027	
150 mm	40 x 9 mm	3,0 mm	BC3028	
185 mm	50 x 9 mm	3,0 mm	BC3029	
nickel 99,5%				
210 mm	55 x 9 mm	3,0 mm	BC3030	
100 mm	35 x 4 mm	2,0 mm	BC3046	
130 mm	40 x 4 mm	2,0 mm	BC3047	
185 mm	50 x 5 mm	2,0 mm	BC3053	
210 mm	55 x 6 mm	2,0 mm	BC3057	
inox AISI 304 18/10 recouvert de Téflon®				
150 mm	40 x 3 mm	1,5 mm	BC3710	
185 mm	50 x 4 mm	2,0 mm	BC3711	
185 mm	50 x 5 mm	2,0 mm	BC3712	
210 mm	55 x 6 mm	2,0 mm	BC3713	

Spatules / lame conique



longueur	lame	Ø fil	référence	Prix HT
inox AISI 304 18/10				
210 mm	50 x 7 mm	2,5 mm	BC3170	
160 mm	50 x 9 mm	3,5 mm	BC3175	
inox AISI 304 18/10 recouverte de Téflon®				
210 mm	50 x 7 mm	2,5 mm	BC3730	

Micro-spatules à manche



• manche en plastique Ø 6 mm

longueur	lame	Ø fil	référence	Prix HT
inox AISI 304 18/10,				
160 mm	40 x 4 mm	2,0 mm	BC3231	
160 mm	40 x 5 mm	2,0 mm	BC3232	

Doubles spatules standard



longueur	lame	Ø fil	référence	Prix HT
acier inox AISI 304 18/10				
130 mm	40 x 9 mm	3,0 mm	BC3100	
150 mm	45 x 9 mm	3,5 mm	BC3101	
185 mm	50 x 9 mm	3,5 mm	BC3102	
210 mm	60 x 11 mm	4,0 mm	BC3103	
250 mm	70 x 11 mm	4,0 mm	BC3104	
300 mm	75 x 16 mm	6,0 mm	BC3105	
350 mm	85 x 16 mm	6,0 mm	BC3106	
400 mm	100 x 16 mm	6,0 mm	BC3107	
500 mm	100 x 20 mm	8,0 mm	BC3108	
nickel pur 99,5%				
130 mm	40 x 9 mm	3,0 mm	BC3110	
150 mm	45 x 9 mm	3,5 mm	BC3111	
185 mm	50 x 9 mm	3,5 mm	BC3112	
210 mm	55 x 11 mm	4,0 mm	BC3113	
250 mm	60 x 11 mm	4,0 mm	BC3114	
acier inox AISI 304 18/10 recouvert de Téflon®				
150 mm	45 x 9 mm	3,5 mm	BC3700	
185 mm	50 x 9 mm	3,5 mm	BC3701	
210 mm	55 x 11 mm	4,0 mm	BC3702	

Spatules à lame flexible



longueur	lame	Ø fil	référence	Prix HT
acier inox AISI 304 18/10				
180 mm	48 x 10 mm	4,0 mm	BC3622	
210 mm	48 x 10 mm	4,0 mm	BC3623	
250 mm	50 x 10 mm	4,0 mm	BC3624	

Spatules plates en nickel pur



longueur	lame	Ø fil	référence	Prix HT
nickel pur 99,5%				
90 mm	20 x 15 mm	1,0 mm	BC3279	
150 mm	32 x 18 mm	1,0 mm	BC3281	

Doubles spatules de Chattaway



longueur	lame	Ø fil	référence	Prix HT
inox AISI 304 18/10				
125 mm	40 x 7 mm	3,0 mm	BC3141	
150 mm	40 x 7 mm	3,5 mm	BC3142	
180 mm	50 x 8 mm	3,5 mm	BC3143	
200 mm	60 x 9 mm	4,0 mm	BC3144	
235 mm	70 x 10 mm	4,0 mm	BC3145	